STATIM[®] 7000

Quick Reference Guide


Operator's Manual should be consulted before using this autoclave

**1 800 667-7733**

**1 800 572-1211**

**+49 (0) 7561 98343-0**



**+1 416 446-4500**

GENERAL USAGE

1


Opening / Closing the Cassette

With thumbs on the thumb pads, push down, pulling up with forefingers until lid opens.


2

Proper Loading


If you notice water beading on the cassette lid, apply Stat-Dri[™] to all internal metal components of the cassette. If the drying performance does not improve, increase the drying cycle time via the User Menu.

NOTE: Always ensure the perforated rack is inserted into the cassette before usage.


3

Inserting the Cassette

Hold the cassette handle in one hand and the carry handle in the other. Place the end of the cassette into the unit and drop the carry handle into its closed position.

Gently push cassette inward until you hear a click. Never force the cassette into the Statim – this could damage the interior components.

CYCLE SELECTION

1

Unwrapped Cycle

This cycle is used to sterilize up to 3.4kg of solid or hollow metal instruments such as pliers, scalers, forceps, dental and phaco handpieces. To select the Unwrapped Cycle, press the Unwrapped Cycle button, then press the START button.

 → 

2

Wrapped Cycle

This cycle is used to sterilize up to 3.4kg of solid or hollow metal instruments which have been sealed in paper or paper and plastic autoclave bags. To select the Wrapped Cycle, press the Wrapped Cycle button, then press the START button.

 → 

3

Rubber and Plastics Cycle

This cycle is used to sterilize up to 0.8kg of unwrapped instruments. To select the Rubber and Plastics Cycle, press the Rubber and Plastics Cycle button, then press the START button.

 → 

4


Air Dry Only Cycle

This is not a sterilization cycle and should not be used if the cassette has been removed from the autoclave. The drying cycle utilizes hot air & has been preset to run for 12 min. To select the air dry only cycle, press the air dry button, then press the START button.

 → 


MAINTENANCE

Condenser Bottle / Waste Bottle

Bottles should be placed below the STATIM 7000.

Replacing the Water Filter

When installed, the water filter has a maximum life of 60 days, depending on the quality of the water. The unit will notify you prior to the expiry or exhaustion of the filter.

Accessing the User Menu

The User Menu can be accessed by using one of two methods:

1. Powering up the unit while pressing the  button or

2. Pressing the  and  buttons simultaneously while the unit displays the (SELECT A CYCLE) message

 = Next option

 = Previous option

 = Select The Option

 = Save & Exit or Return To The Main Menu

MAINTENANCE SCHEDULE

To ensure trouble-free performance, operators must follow this preventative maintenance schedule.

DAILY		WEEKLY		EVERY 6 MONTHS																	
Run Air Dry Cycle	It is recommended to run an Air Dry Cycle at the end of every work day to allow any residual moisture remaining in the system to dry.	Cassette	Wash the interior of the cassette with dishwashing soap or a mild detergent that does not contain chlorine.	Biological and / or Air Filter	Replace every 500 cycles or six months (whichever is first)																
Water Reservoir	Check the reservoir for dirt. Clean and rinse if necessary.		Scrub the inside with a cleaning pad designed for use with Teflon™ – coated surfaces.																		
	For ophthalmic use, drain at the end of every workday, leave empty, and refill at the start of the next workday.		After removing all traces of the detergent, treat interior surfaces of the cassette with the STAT-DRI drying agent to enhance the drying process. Order more STAT-DRI from SciCan quoting 2OZPLUS, 8OZPLUS, or 32OZPLUS.																		
Waste Bottle	Empty the waste bottle every time you refill the water reservoir.	Biological and / or Air Filter	Check the filter for dirt and moisture. Replace if dirty. Call for service if wet.																		
	Check to ensure fittings are secure	Condenser bottle	The condenser bottle should be cleaned and refilled to the maximum.																		
Every 250 cycles	Lubricate cassette and reservoir O-Rings Seals or lubricate when cassette is difficult to engage																				
EVERY YEAR																					
Cassette Seal / O-Rings		Replace every 1,100 cycles or 12 months (whichever is first).																			
SPARE PARTS																					
		<table><tr><td>• Water Filter (6 pack)</td><td>SCWF6</td></tr><tr><td>• Water Filter (single)</td><td>SCWF1</td></tr><tr><td>• Cassette Seal Kit</td><td>01-110295S</td></tr><tr><td>• Mat – Silicone 7000</td><td>01-110754S</td></tr><tr><td>• Bottle Overflow To Drain</td><td>01-110300S</td></tr><tr><td>• STAT-dri 8 oz.</td><td>8OZPLUST</td></tr><tr><td>• STAT-dri 32 oz.</td><td>32OZPLUS</td></tr><tr><td>• By-pass Water Filter</td><td>SCWFBP</td></tr></table>				• Water Filter (6 pack)	SCWF6	• Water Filter (single)	SCWF1	• Cassette Seal Kit	01-110295S	• Mat – Silicone 7000	01-110754S	• Bottle Overflow To Drain	01-110300S	• STAT-dri 8 oz.	8OZPLUST	• STAT-dri 32 oz.	32OZPLUS	• By-pass Water Filter	SCWFBP
• Water Filter (6 pack)	SCWF6																				
• Water Filter (single)	SCWF1																				
• Cassette Seal Kit	01-110295S																				
• Mat – Silicone 7000	01-110754S																				
• Bottle Overflow To Drain	01-110300S																				
• STAT-dri 8 oz.	8OZPLUST																				
• STAT-dri 32 oz.	32OZPLUS																				
• By-pass Water Filter	SCWFBP																				

 Best Practices: SciCan recommends the final user carefully choose the most appropriate sterilization cycle according to the recommendations of their leading infection control authorities and local regulatory guidelines / recommendations.

Your Infection Control Specialist[™]

www.scican.com

STATIM is a registered trademark and Stat-Dri, Dri-Tec[®] and “Your Infection Control specialist” are trademarks of SciCan Ltd.

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