M3 UltraFast® Automatic Sterilizer

[English](#)

[Español](#)

[Français](#)

For Models:
M3 (-001 thru -004)

User Guide


Product Information

Dealer:

Date of Purchase:

Model / Serial Number:

Midmark Authorized Service Company:

SA1794i

Table of Contents

Important Information

Component Location	4
Safety Instructions	5
Intended Use	5
Electromagnetic Interference	5
Safety Symbols	6
Transportation / Storage Conditions	6
Operating Environment	7
Electrical Ratings / Requirements	7
Sterilization Monitoring Guidelines	8

Installation

Location Requirements	9
Power Cord Connection	10
External Condensing Tank Connections	11
External Condensing Tank Draining Procedures	13

Operation

Quick Reference	14
Power Switch	15
Filling the Reservoir	15
Qualification Testing	16
Loading the Tray	16
Cycle Selection	20
Post Sterilization Processing	22
'Additional Heat' Cycle	22
Adjusting the Dry Time	23

Accessories

Accessory Chart	23
-----------------------	----

Maintenance

Maintenance Messages	23
Daily Care	24
Periodic Maintenance	25
Monthly Maintenance	27
Extended Use Maintenance	28

Printer (optional)

Printer Harness Connection	29
Printer Readout	30
Paper Roll Installation	31
Printer Ribbon Cartridge Replacement	32

Troubleshooting

Error Codes	33
User Diagnostic Mode	34
Calling for Service	35

Specifications


Specification Information	36
Water Purity Specifications	37

Warranty Information

Scope of Warranty	38
-------------------------	----

Important Information

Component Location

SA1795i

Safety Instructions

The primary concern of Midmark is that this equipment is operated and maintained with the safety of the patient and staff in mind. To assure safe and reliable operation:

- Read and understand this manual before attempting to install or operate the sterilizer.
- Assure that the appropriate personnel are informed on the contents of this manual.
(This is the responsibility of the purchaser).
- Assure that this manual is located near the sterilizer, or if possible, permanently affixed to the sterilizer.

Intended Use

The M3 UltraFast® Automatic Sterilizer can be used in medical, dental, and veterinary offices, hospitals, clinics, nursing homes, laboratories, and other facilities to sterilize heat and moisture stable reusable items (including dental handpieces) that are compatible with steam sterilization.

Refer to 'Loading the Tray' & 'Cycle Parameters' in this manual for detailed information.

Electromagnetic Interference

The Midmark M3 is designed and built to minimize electromagnetic interference with other devices.

However, if interference is noticed between another device and this sterilizer:

- Remove interfering device from room
- Plug sterilizer into a dedicated circuit
- Increase separation between sterilizer and interfering device
- Contact Midmark if interference persists

Safety Symbols

DANGER
Indicates an imminently hazardous situation which will result in serious or fatal injury. This symbol is used only in the most extreme conditions.

WARNING
Indicates a potentially hazardous situation which could result in serious injury.

Caution
Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices

Equipment Alert
Indicates a potentially hazardous situation which could result in equipment damage.

Note
Amplifies a procedure, practice, or condition.


Shipping Symbols

Caution Shipping Damage

Keep dry

Proper shipping orientation

Fragile

2 Person Lift

Transportation / Storage Conditions

Equipment Alert
All water must be removed from the reservoir before transporting or storing at +32°F (0°C) or below.

Ambient Temperature Range:..... -22°F to +140°F (-30°C to +60°C)
Relative Humidity:..... 10% to 90% (non-condensing)
Atmospheric Pressure: 49.6 kPa to 106.4 kPa (7.2 psi to 15.4 psi)

Operating Environment

Equipment Alert

*Unit should be allowed to reach room temperature before operating.
Failure to do so could result in damage.*

Ambient Temperature Range:..... +68°F to 104°F (+20°C to 40°C)

Relative Humidity:..... less than 80% (non-condensing)

(Pollution Degree 2, in accordance to IEC664)

Normal Operating Altitude: less than 9842 ft. (3000 m) above sea level

- Approved for indoor use only
- Environment should be relatively dust-free

Electrical Ratings / Requirements

Note

To ensure unit is properly grounded, it must be connected to a matching grounded, dedicated, correctly polarized receptacle.

WARNING

Use 104-127 VAC, 50/60 HZ alternating current only for 115 VAC rated models and 207-253, 50/60 HZ alternating current only for 230 VAC rated models. Failure to do so could result in electrical shock to personnel and will result in damage to sterilizer.

M3 (115V model):115 VAC, 50/60 Hz, 12 amp

Max. Power Consumption:1400 Watts

Requires*:.....Dedicated supply circuit rated at 120 VAC, 50/60 Hz, 12 amp

M3 (230V model):.....230 VAC, 50/60 Hz, 6 amp

Max. Power Consumption:1400 Watts

Requires*:.....Dedicated supply circuit rated at 230 VAC, 50/60 Hz, 6 amp

**Power source must have over voltage limits less than 1500 watts from mains to ground.*

(Installation Category II in accordance to IEC 664)

Sterilization Monitoring Guidelines

Note

The information below is provided for reference only. Contact appropriate state / local agencies for specific sterilization guidelines for your office. Additional information on infection control is available from the Centers for Disease Control and Prevention (CDC), Organization for Safety and Asepsis Procedures (OSAP), and the American Dental Association (ADA).

Physical Monitors

Temperature and pressure measuring devices can help detect sterilizer malfunctions. The sterilizer's control system aborts the cycle and displays a message if physical conditions go outside established limits. The optional printer can be used to create a record of each load's actual cycle time, temperature, and pressure.

Note

Use only FDA cleared chemical & biological indicators designed for steam sterilization that are compatible with the particular sterilization cycle temperature and exposure time being monitored. Use sterility monitors with each sterilization load. If a sterilization cycle is terminated prematurely, reprocess instruments to ensure sterility of the load. Process the load according to your regular practice, placing indicators near the handle side of tray. Follow manufacturer's instructions for proper disposal of used indicators.

Chemical Indicators

Chemical indicators are designed to verify that conditions in the sterilizer chamber were adequate to achieve sterilization. They do not validate that a processed item is sterile. If a chemical indicator shows a failure, items in that load are considered non-sterile. Potential causes for sterilization failures include: improper cleaning, packing, loading, or a sterilizer malfunction. Determine the cause of any sterilization failure, and remedy the situation before running the next cycle. Only FDA cleared chemical indicators labeled for use with the nontraditional steam sterilization cycle parameters, e.g. temperature and exposure time, of the M3 Sterilizer should be used for monitoring the three M3 UltraFast® cycles. Follow the chemical indicator's instructions for proper storage, use, interpretation, and disposal.


Biological Indicators

Biological indicators are microbiological devices designed to accompany items being sterilized to monitor adequacy of the sterilization process. If a biological indicator shows a failure, items in that load are considered non-sterile. Potential causes for sterilization failures include: improper cleaning, packing, loading, or a sterilizer malfunction. Determine the cause of any sterilization failure, and remedy the situation before running the next cycle. Only FDA cleared biological indicators labeled for use with the nontraditional steam sterilization cycle parameters, e.g. temperature and exposure time, of the M3 Sterilizer should be used for monitoring the three M3 UltraFast® cycles. Follow the biological indicator's instructions for proper storage, use, interpretation, and disposal.

Installation

Location Requirements

(Allow clearance on both sides)

SA1796i

Support Surface

- Material should be water-resistant material
(Ex. laminate, stainless steel, stone, etc.)
- Surface must be level to ensure proper operation.
- Surface should meet minimum dimensions listed below:

Dimensions

Depth (front to back).....24 in. (61 cm)

Width (side to side).....22 in. (56 cm)

Clearance Requirements

To ensure proper air circulation, and to allow access to the reservoir fill port and drain coupling, adhere to the minimum clearance requirements listed below.

Clearance Requirements

Back of Unit - Back Wall	4 in. (10 cm)
Front Sterilizer Feet - Front of Support Surface	4 in. (10 cm)
Side of Unit - Side Wall	2 in. (5 cm) each side
Distance Above Unit*	2 in (5 cm)*

* The minimum clearance for proper air circulation is listed.
However, be sure to allow access to the reservoir fill port located on top of the sterilizer.

Power Cord Connection

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture with oxygen, air, or nitrous oxide.

Clarification: Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.

WARNING

Check the serial number label on back panel of sterilizer to verify voltage rating for the unit. Failure to connect sterilizer to an appropriate power supply could result in damage to the unit, and electrical shock to personnel.


Equipment Alert

For optimal performance, allow sterilizer to reach room temperature before operating.

To connect the power cord...

- A) Plug power cord into receptacle on back of sterilizer.
- B) Plug power cord into a properly polarized and grounded receptacle rated for a minimum of 15 amps. A dedicated circuit only used for the sterilizer is recommended.

SA1816i

External Condensing Tank Connections


Note

External Condensing Tank is not needed if M3 UltraFast® is connected to the VistaCool direct-to-drain thermal reduction system.

Equipment Alert

The tank hose must be properly connected (and not kinked). Improper connection will cause water / steam leaks or a sterilizer malfunction if water flow to the tank is restricted.


To connect the condensing tank...

- If connected, disconnect sterilizer power cord.
- Push tank hose into compression fitting on lid. Tighten nut. (If step B is difficult - see NOTE)
- Connect tank hose to back of sterilizer.
- Fill condensing tank with tap water to minimum water level indicator line.
- Secure the lid on the condensing tank by aligning the tabs and rotate lid clockwise to secure.

Note

If tubing is difficult to install...

- Remove compression nut & two ferrules.
- Install nut and ferrules onto tank hose as shown. (Large & small ferrules must be positioned as shown)
- Insert hose into tank fitting, then tighten nut.


Compression Fitting

Min. Water Level Indicator Line

SA1826i

External Condensing Tank Connections - continued

To connect the condensing tank - continued...

- F) Place tank on a level surface, preferably below the sterilizer but in no case should it be higher than the sterilizer support surface.
- G) Connect sensor plug.
- H) Connect sterilizer power cord.

Note

Clearance...

Maintain a minimum of 6 inch clearance above the condensing tank for proper steam ventilation. If enclosed in a cabinet the support surface and surrounding surfaces, should be protected with a water resistant material (e.g. plastic, laminate, stainless steel, etc.). If enclosed in a cabinet, it is recommended that the door be vented to avoid heat, moisture build up and potential damage to the inside of cabinet.

SA1827i

External Condensing Tank Draining Procedure


Caution

Water that is discharged to external condensing tank can be **VERY HOT**; person emptying pitcher should allow the temperature to cool. Always use carrying handle and use caution when emptying.


To remove lid...

- A) Rotate lid counterclockwise.
- B) Lift lid assembly out of condensing tank pitcher.

To drain pitcher...

- A) Grasp pitcher by handle.
- B) Pour cool water into an approved drain.

To replace lid...

Note: Lid can be reinstalled in one of four locations to facilitate safe water disposal.

- A) Insert lid assembly into pitcher.
- B) Align the four tabs for correct positioning.
- C) Rotate lid clockwise to secure.

SA1828i



Operation

Quick Reference

(Detailed instructions for each step are outlined in the following pages of the Operation section).

Basic Operation...

- A) Turn power switch ON (I).
- B) Fill reservoir.
- C) Load tray.
- D) Press desired cycle button.
- E) Press <Start> button.

Note

There is a 10 minute (dry) heat cycle that can be used to pre-heat the chamber, or for additional drying time. Refer to: '**Additional Heat**' Cycle in this section.


SA1799i

Unwrapped Pouches Low Temp

Cycle Buttons

Power Switch

The power cord must be connected and the power switch must be ON (I) for the sterilizer to operate.



SA1800I

Filling the Reservoir

To fill reservoir...

Pour one (1) gallon of distilled water into fill port.
Do not fill above lower lip of fill port.

Equipment Alert

Use distilled water or water that meets the referenced water purity specifications. Failure to comply may result in sterilizer malfunction and/or premature failure due to excessive corrosion.

SA1801I


Qualification Testing

Your sterilizer should be tested after sterilizer installation, malfunctions, relocation, major repairs, and after sterilization process failure. Qualification testing should be performed prior to placing the sterilizer in service. If multiple cycles types are used, e.g. “Pouches” and “Low Temp” each cycle type should be qualified. Qualification testing should include at least one Biological Indicator (BI) (sometimes referred to as Spore Tests) and one Chemical Indicator (CI). The test pack should be performed with items routinely processed and considered to be the most difficult to sterilize. Additional items should be placed in the chamber along with the Biological Indicator and Chemical Indicator so that chamber is fully loaded (don't exceed the maximum capacities listed in the tables under “Guidelines for Loading” in this manual). Three consecutive test runs, for each cycle type tested, with negative results from the BIs, and the appropriate readings from all physical monitors and chemical indicators demonstrating complete sterilization, provide verification that the sterilizer has been properly installed (or reinstalled after relocation) or repaired to the manufacturer's specifications and that it will function effectively in the facility in which it is installed. All items processed during qualification testing should be quarantined until the results of the biological testing for all three test runs are available.

Loading the Tray

Types of Items

Before placing any instrument in the M3 UltraFast®, check with the instrument manufacturer to be sure the materials are compatible with steam sterilization, and to verify the acceptability of sterilization parameters.

The M3 is designed to sterilize the following:

- High & low speed handpieces
- Metal instruments
- Rubber / plastic devices (ex. suction cannulas, impression trays, etc.)
- Wrapping / bundling materials (ex. CSR wrap, instrument pouches, etc.)
- Cassettes (Hu-Friedy Signa-Stat [6.5" x 10.5" x 1.25"] or smaller)
- Surgical instruments (ex. ophthalmologic instruments)

Equipment Alert

Do not sterilize items composed of any of the following materials in the M3:

- Corrosion sensitive metal (ex. carbon steel, iron, etc.)
- Fragile items susceptible to breaking under pressure / high temperature
- Liquids
- Biomedical waste
- Textiles (including towels, gauze, etc.)
- Plastics that may break down or produce residue when exposed to steam / high temperatures.

Examples

Polyethylene	Styrene	Cellulosics	ABS
PVC	Textiles	Acrylic (Plexiglass™)	
PPO (Noryl™)	Latex	Neoprene	

Loading the Tray - continued

Immediate Use Sterilization


The M3 is capable of immediate use sterilization - sterilizing unwrapped instruments for immediate use. Please consider the following when choosing whether or not to sterilize your instruments unwrapped:

- The sterility of unwrapped instruments is compromised upon exposure to a non-sterile environment. Follow CDC guidelines for using unwrapped, sterilized instruments.
- Due to the sensitive nature of some types of surgery (including, but not limited to ophthalmological), instruments used in such procedures must be wrapped or pouched in order to reduce their exposure to sterilization process residues. The water reservoir should also be drained and refilled with fresh distilled water on a daily basis when processing instruments for these procedures on a routine basis.

Pouching and Wrapping Items

The M3 is capable of sterilizing pouched or wrapped items.

- When pouching or wrapping items, use only sterilizer pouches and wraps that have been cleared by the FDA and labeled for use with the nontraditional steam sterilization cycle parameters, e.g. temperature and exposure time, of the M3 sterilizer. Follow the manufacturer's instructions for use.
 - When using cassettes in the M3 follow the manufacturer's instructions for use.
 - Pouched items to be sterilized should be placed lengthwise with plastic side of pouch facing up in the M3 UltraFast® tray.
- The pouches may overlap slightly, but items must not be layered. Refer to diagram below.

SA1803i

Caution

Failure to comply with these guideline may adversely affect sterilization and/or drying.

Loading the Tray - continued

Load Size

The M3 UltraFast® can accommodate loads weighing up to **2.4 lbs (1.1 kg)**.

[Note: This is the weight of the contents in the tray (ex. instruments, cassettes, pouches, etc.).

The weight of the tray itself has already been accounted for].

WARNING

Do not overload the sterilizer tray. Failure to allow adequate space around items for circulation will compromise sterilization and drying.

Use the table below as a general guideline for weights of commonly used items. Consult manufacturer's specifications for the exact weight of any particular instrument.

Item Description	Weight*	
	lbs.	kg
Scissors	0.066	0.030
Dental Scalars	0.044	0.020
Forceps	0.033	0.015
Dental Handpiece	0.121	0.055
Suction Cannula	0.022	0.010
Plastic Mouth Mirror	0.018	0.008
Impression Tray	0.033	0.015
Plastic X-Ray Positioning Ring	0.044	0.020
Hu-Friedy Signa-Stat Cassette	1.500	0.680

Packing the Tray

(*actual weights may vary)

Equipment Alert

Do not use towels or packaging which contains chlorine bleach residue. Chamber and/or tray rusting or discoloration may occur. The life of the sterilizer may be shortened significantly.

Caution

Clean and dry instruments thoroughly before placing them into tray. Improper cleaning may result in non-sterile instruments or damage to the unit. Follow instrument manufacturer's guidelines and CDC recommendations for handling and cleaning instruments prior to sterilization.

In addition to total load weight outlined above, all items must be processed in accordance with Centers for Disease Control and Prevention (CDC), 'Guidelines for Infection Control in Dental Healthcare Settings' - 2003, MMWR 2003; 52 (no. RR-17), which states:

"Items to be sterilized should be arranged to permit free circulation of the sterilizing agent (e.g., steam, chemical vapor, or dry heat); manufacturer's instructions for loading the sterilizer should be followed."


- All items must fit in M3 UltraFast® tray.
- Loaded tray must slide into chamber opening without scraping.
- Items must not touch one another.
- Pouched items should be loosely packed.
- Pouches may overlap slightly, but items must not be layered.

Caution


Failure to comply with these guidelines may adversely affect sterilization and/or drying.

Loading the Tray - continued

Cycle Selection

(The parameters for each cycle are outlined on the following page.)

Cycle Selection - continued

Cycle Parameters

(Before sterilizing any items in the M3, refer to **Loading the Tray** in this section.)

CYCLE	STERILIZATION PARAMETERS	DRY TIME	ITEMS TO BE STERILIZED
 Unwrapped	Temperature: 270°F (132°C) Pressure: 27.1 psi (186 kPa) Time: 3:30 Minutes	Time: 25 Minutes	- Dental instruments / handpieces loose on a tray. - Items manufacturers recommend for exposure at 270°F (132°C), loose on tray.
 Pouches	Temperature: 270°F (132°C) Pressure: 27.1 psi (186 kPa) Time: 5:30 Minutes	Time: 30 Minutes	- Dental instruments / handpieces in pouches, wrapped, or in a wrapped cassette. - Items manufacturers recommend for exposure at 270°F (132°C), in pouches, wrapped, or in a wrapped cassette.
 Low Temp	Temperature: 250°F (121°C) Pressure: 15.0 psi (104 kPa) Time: 20:00 Minutes	Time: 50 Minutes	- Rubber or plastic dental devices, dental instruments / handpieces loose on a tray, in pouches, wrapped, or in a wrapped or unwrapped cassette. - Items manufacturers recommend for exposure at 250°F (121°C), loose on tray, in pouches, wrapped, or in a wrapped or unwrapped cassette.

Post-Sterilization Processing


After sterilization is complete, all items must be handled in accordance with accepted and documented standards, such as the Centers for Disease Control and Prevention (CDC) document, 'Guidelines for Infection Control in Dental Healthcare Settings' - 2003, MMWR 2003; 52 (no. RR-17), as well as any local requirements that may apply.

Qualified personnel responsible for infection control should prepare a protocol for handling sterilized items. This protocol should be followed by all personnel responsible for handling sterilized items.

'Additional Heat' Cycle

The Additional Heat Cycle activates the dry heaters for ten minutes.

This cycle can be used to pre-heat the chamber at the beginning of the workday, or for extended drying time at the end of a cycle.

SA1804i

To pre-heat chamber prior to running a cycle...


- A) Press <Start> button when 'SELECT CYCLE' appears on display.
- B) During the ten minute pre-heat mode, 'ADDITIONAL HEAT' will flash on the display.
- C) When 'ADDITIONAL HEAT' stops flashing, press desired cycle button, then press <Start>.

For extended drying time at the end of a cycle...

- A) Press <Start> button when 'SELECT CYCLE' appears on display.
- B) During the ten minute drying mode, 'ADDITIONAL HEAT' will flash on the display.

Adjusting the Drying Time

The M3 allows the operator to adjust the drying time from 20 - 60 minutes using 1 minute increments for the three pre-programmed cycles.

SA1805i

To adjust the drying time for a pre-programmed cycle...

A) After pressing desired cycle button, press <P> button.

[Display will show current setting. (ex. DRY TIME: 30 MINUTES)]

B) Press the <+> / <-> buttons to increase / decrease the drying time.

C) Press the <P> button to save your changes.

(Pressing the <Stop> button cancels the changes & returns to last saved setting.)

Accessories

Accessory Chart		
Description	Part Number	Intended Use
Printer	9A401001	Prints sterilizer cycle information and data.
Door Tray	9A402001	Used to stage cycles.
Top Cover Protector	9A404001	Used to protect painted top cover from damage.

Maintenance


Maintenance Messages

To assure correct operation and maximum sterilizer life, the M3 provides the operator with reminders when it's time to perform operator maintenance. After the M3 is powered ON for 7, 14, and 21 days, a message "Perform Periodic Maintenance" will be displayed. After 28 days, a "Perform Monthly Maintenance" message will be displayed. Refer to the appropriate maintenance instructions in this manual. The maintenance reminders are removed from the display when a cycle is started. If power is turned OFF, the timer will reset, initiating a new cycle of messages.

Daily Care


- **Clean External Surfaces / Tray & Chamber**

- A. Wash unit according to your facility's procedure for clinical contact surfaces noting the following:
(Use only quaternary disinfectants to disinfect unit. Staining, pitting, discoloration, or softening could occur if phenolic, iodophor, or glutaraldehyde-based disinfectant is used on plastic surfaces of the unit. Also, use of alcohol or aerosol spray cleaner / disinfectant containing substantial amounts of alcohol in the formula can damage the faceplate.)
- B. Wring excess solution from cloth.
- C. Using soft cloth, wipe all external surfaces.
- D. Follow the instructions provided with the cleaner / disinfectant used regarding rinsing and drying of external surfaces.

WARNING

Check to assure the chamber filter is not blocked by debris and that it is properly seated; in the bottom hole in the back of the chamber. Failure to make sure the filter is in place and clear of debris could result in serious personal injury and or equipment damage.

A circular inset showing a close-up of the chamber filter, which is a cylindrical mesh screen.


Chamber Filter

- **Clean Chamber Filter**

Check that the chamber filter is free of debris and properly positioned in the hole at the back of the sterilizer chamber. If the filter is blocked by debris follow the procedures under "Monthly Maintenance" to remove and clean the filter.


- **Clean Door Gasket / Mating Surface**

- A. Wash with a damp cloth.
- B. Inspect gasket for damage.
- C. Replace gasket if necessary.

SA1806j

Periodic Maintenance

- **Drain / Refill Reservoir (As Needed)**

- Remove drain hose from storage location.
- Place open-end of drain hose into container or sink.
- Connect adapter-end of drain hose to coupling as shown.
- Once water has drained, press release lever and remove hose.
- Return drain hose to storage location.
- Refill reservoir with distilled water.

Equipment Alert

Use distilled water or water that meets the referenced water purity specifications. Failure to comply may result in sterilizer malfunction and/or premature failure due to excessive corrosion.

Periodic Maintenance - continued

• Empty / Clean External Condensing Tank (As Needed)

Note

This process is not necessary if the sterilizer is connected to a direct-to-drain thermal reduction system.

- A. Empty water from tank. (Do not reuse water!)
- B. Clean tank with diluted bleach solution (1/4 cup bleach : 1 gallon water) and a brush.
- C. Rinse tank thoroughly.
- D. Refill tank to minimum water level indicator line.

SA1808i

Monthly Maintenance

- **Remove & Clean Filter**


Wash with mild soap solution to remove debris. Rinse with distilled water.
(Use a stiff brush to scrub, or place in ultrasonic cleaner if necessary.)

- **Clean Condensing Tank Level Sensors**

Note


This process is not necessary if the sterilizer is connected to a direct-to-drain thermal reduction system.

Clean two sensors with mild soap solution, then wipe dry.


SA1809i

Monthly Maintenance - continued

To remove filter...

- A) Slide notch in filter tool (provided) over filter.
- B) Pull out to remove.

To install filter...

- A) Position filter in notch of tool as shown.
- B) Align filter with hole in back of chamber*, and press in gently.
- C) Flip tool over, then gently press filter in to secure.

*Tip: Slide tool along bottom of chamber to align filter with hole.

SA1810i

Extended Use Maintenance

The M3 is designed and tested to provide exceptional reliability throughout its service life. However, like all electro-mechanical devices it is subject to wear and degradation with use.


To ensure the integrity, performance, and safety of all major components it is the responsibility of the user to have the sterilizer performance / operation verified by a Midmark Authorized Service Provider after 8 years or 20,000 cycles of use, whichever comes first. After 8 years or 20,000 cycles of use an annual inspection by a Midmark Authorized Service Provider is recommended.

Printer (optional)

Printer Harness Connection

To connect printer harness...

- A) Turn sterilizer power switch OFF (O).
- B) Connect printer harness as shown.
- C) Turn sterilizer power switch ON (I).

SA1813I

Printer Readout

Selected Cycle

Total # of cycles
run on unit

During Sterilization phase,
the printer records chamber
temperature & pressure in
30-second increments.

Indicates Venting
phase initiated

Set Temperature & Time

Date

During Heat-Up phase, the printer
records chamber temperature &
pressure in 30-second increments.

Summary of Sterilization phase

Indicates Drying phase initiated
(Printer records drying time
in 5 minute increments).

Paper
Feed Button

If printer readout becomes light ...

- Separate lid / printer from housing.
- Press Paper Feed button for three seconds.
(This advances the ribbon inside the cartridge).
- Repeat step B as necessary*.

*Note: If pressing the Paper Feed button fails to correct
the problem, replace the printer cartridge.


```

BEGIN
UNWRAPPED
270 deg F 3.5 MINUTES

MM - DD - YYYY
__ - __ - __
HH : MM
__ : __
TOTAL CYCLES 140

HEATING

MM:SS  degF  PSI
0:00  151.5  0.1
0:30  193.1  0.2
1:00  218.1  3.8
1:30  243.2  11.9
2:00  244.6  13.5
2:30  268.5  27.8

STERILIZING

MM:SS  degF  PSI
0:00  271.0  28.8
0:30  274.8  31.6
1:00  272.9  29.2
1:30  275.1  25.8
2:00  273.7  30.8
2:30  273.8  30.0
3:00  273.5  29.3
3:30  275.2  25.8

MIN TEMP  271.0  degF
MAX TEMP  275.7  degF
MIN PRESS  23.2  PSI
MAX PRESS  32.1  PSI

VENTING CHAMBER

DRYING

MM:SS
0:00

TOTAL PROCESS TIME
9:50

COMPLETE
    
```


SA1814i

Paper Roll Installation

To install paper roll...




- A) Loosen thumbscrew, then separate lid / printer from housing.
- B) Slide paper into slot as shown*.
- C) Press paper feed button until approx. 2 in. (5 cm) of paper is fed thru printer.
- D) Pull paper thru slot in lid.
- E) Position paper roll / spindle in slots.
- F) Reassemble lid / printer / housing.

Printer Ribbon Cartridge Replacement

To remove old printer ribbon cartridge...

- A) Turn sterilizer power switch OFF (O).
- B) Separate lid from printer.
- C) Press down on side of printer cartridge labeled "EJECT".

To install new printer ribbon cartridge...

- A) Install new cartridge as shown* (it will "snap" into place).

*Note: The printer paper must be between the cartridge ribbon & the metal plate.

- B) Turn Adjustment Knob clockwise until ribbon is tight.

Troubleshooting

Error Codes

If a malfunction is detected during a cycle, a numeric error code will appear on the display panel. Use the chart below to diagnose and correct the most common, maintenance-related error codes. If you encounter an error code not identified below, follow the instructions on the display panel. If error code persists, contact your authorized service provider.

Example:

Error Code	Probable Cause	Corrective Action
C010	Sterilizer lost power during cycle.	Press <STOP> button, then restart cycle.
C100-series (all C101, C102, etc.)	<STOP> button was pressed during the cycle.	Press <STOP> button, then restart cycle.
C231, C232	Not enough water in reservoir to complete the cycle.	Fill reservoir with distilled water, or water that meets the referenced water purity specifications.
C441, C442	External condensing tank is full.	Empty external condensing tank.
C533, C633	Water pump needs primed.	Put the sterilizer in the User Diagnostic Mode and prime the pump. (see following page)

Troubleshooting - continued

User Diagnostic Mode

The User Diagnostic Mode is used:

- To set "English" or "Metric" units on the display panel,
- To retrieve the last five (5) error codes stored in the unit memory.
- To prime the water pump if needed.

To activate User Diagnostic Mode...

- A) Turn Power Switch OFF (O).
- B) Press and hold the <START> button.
- C) Turn Power Switch ON (I).
- D) Press the <START> button when the display shows "USER DIAGNOSTIC".

To change Display Units...

- A) Put the unit in User Diagnostic Mode.
- B) Press the <P> button to select units.
- C) Press the <+> button to change the temperature and pressure display from English-to-Metric or Metric-to-English units. (Factory default setting is English).
- D) Press the <START> button to continue.
- E) Turn the power switch OFF (O) to exit User Diagnostic Mode.

To retrieve the five (5) most recent error codes...

- A) Put the unit in User Diagnostic Mode.
- B) Press the <STOP> button to recall errors.
- C) The last five (5) error codes will be displayed.
- D) Press the <START> button to return to the User Diagnostic Mode display.
- E) Turn the power switch OFF (O) to exit User Diagnostic Mode.

To prime the sterilizer pump...

- A) Put the unit in User Diagnostic Mode.
- B) Press the <START> button to start the pump priming progress.
The unit will automatically cycle through a pre programmed priming cycle...
 - Closing the sterilizer door.
 - Heating the boiler.
 - Cycling the pump ON and OFF until the pump is primed.When finished the 2nd line of the display will show "PRIMING COMPLETE".
- C) Press the <START> button to return to the User Diagnostic Menu.
- D) Turn the power switch OFF (O) to exit User Diagnostic Mode.

Calling for Service

Note

Please mark down any displayed Code(s) and be sure to relay this information to the service technician.

Contact your **Midmark Authorized Dealer**, or log on to www.midmark.com/technical-library.
Model and serial number information will be required when calling for service.

To contact Midmark directly:

1-800-MIDMARK (1-800-643-6275) or 937-526-3662
8 am to 5 pm EST (Monday thru Friday)
[excluding standard U.S. holidays]

Specifications

Fuse Ratings:

<i>115 VAC Unit</i>	
F1	15 Amp, 250 V, Fast Acting, 1/4" x 1 1/4"
F2	0.25 Amp, 250 V, Slo-blo, 1/4" x 1 1/4"
<i>230 VAC Unit</i>	
F1	8 Amp, 250 V, Fast Acting, 5 x 20 mm
F2	0.125 Amp, 250 V, Slo-blo, 5 x 20 mm

Certifications:

ASME Boiler & Pressure Vessel Code, Section VIII, Division 1
Canadian Registration Number Available
UL 61010-1, 2nd Edition
IEC 61010-2-040
CAN/CSA-C22.2 No. 61010-1, 2nd Edition
FCC Part 15, Sub-part B

Physical Dimensions:

Overall Length:	21 in. (53.3 cm)
Overall Width:	17.8 in. (45.2 cm)
Overall Height:	7.1 in. (18 cm)
Shipping Carton Length:	25 in. (63.5 cm)
Shipping Carton Width:	22 in. (55.9 cm)
Shipping Carton Height:	16.6 in. (42.2 cm)
Counter Area:	24 in. (61 cm) deep x 22 in. (55.9 cm) wide
Chamber Volume:	0.49 gal (1.8 liter)

Weight:

Empty Reservoir:	71 lbs. (32.2 kg)
Full Reservoir:	80 lbs. (36.3 kg)
With Shipping Carton:	80 lbs. (36.3 kg)

Water Reservoir Capacity: 1.20 gal (4.5 liter)

Pressure Relief Valve Setting: 40 PSI (275.8 kPa)

Chamber Pressure:

@ 270°F (132°C) 27.1 psi. (186.2 kPa)

Water Purity Specifications

Water Purity Specification Chart	
	AAMI ST-46 (ref.)
Evaporate Residue	≤ 15 mg/l
Silica	≤ 2 mg/l
Iron	≤ 0.2 mg/l
Cadmium	≤ 0.005 mg/l
Lead	≤ 0.05 mg/l
Rest of heavy metals, excluding Iron, Cadmium, & Lead ¹	≤ 0.1 mg/l
Chloride	≤ 3 mg/l
Phosphate	≤ 0.5 mg/l
Conductivity (at 20°C) ²	≤ 50 µs/cm
pH	6.5 to 8
Appearance	Colorless, clean, w/o sediment
Hardness	≤ 0.1 mmol/l (.585 gr/gal)

Warranty Information

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the first retail purchaser of a product and is not transferable or assignable.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be one (1) year for all warranted products and components.

OBTAINING WARRANTY SERVICE

Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com; by phone at 1-800-MIDMARK; by facsimile at 1-800-365-8631; or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this warranty.

EXCLUSIONS

This warranty does not cover, and Midmark shall not be liable, for the following:

- (1) defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) products which are not installed, used, and properly cleaned and maintained as required in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product;
- (3) products considered to be of a consumable nature;
- (4) accessories or parts not manufactured by Midmark;
- (5) charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) costs and expenses of routine maintenance and cleaning; and
- (7) representations and warranties made by any person or entity other than Midmark.

Warranty Information (continued)

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER:

MIDMARK'S ONLY OBLIGATION UNDER THIS WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

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