













# Preventative Maintenance Checklist

Frequency	Preventative Maintenance	Item	Completed by	Date	Order Link
Daily	<p><b>Clean the External Surfaces</b></p> <p>Clean all external sterilizer covers with a slightly damp cloth moistened with water. <i>Notice: Never use detergents or abrasive products.</i></p>				 <p>Optim1 6"x7" Wipes</p>
Monthly or 50 Cycles	<p><b>Clean the Door Gasket &amp; Chamber Face Side</b></p> <p>1 Clean the door gasket and the outer edge of the chamber with a non-abrasive cloth moistened with clean water. 2 Rinse with clean water. <i>Note: when the seal is new it might be necessary to hold the door gently closed at the sterilization start.</i></p>	<p>OPTIM1 CLEANER</p> 			 <p>W&amp;H Lexa Chamber Filters</p>
	<p><b>Clean the Chamber &amp; Chamber Accessories</b></p> <p>1 Remove the trays and the chamber rack. 2 Clean the chamber with a damp sponge and a mild detergent solution paying attention not to bend or damage the temperature probe inside the sterilizer chamber. 3 Rinse with water. 4 Clean the trays and the chamber rack with a damp sponge and a mild detergent solution. 5 Rinse with water. 6 Reposition all pieces of the chamber accessories properly. <i>Note: The trays and the tray holder may also be cleaned in a washer disinfectant.</i></p>	<p>LEXA CHAMBER FILTERS</p> 			 <p>W&amp;H Lexa Door Seal</p>
	<p><b>Clean the Chamber Filters</b></p> <p>1 Allow the sterilization chamber to cool down. 2 Empty the sterilizer chamber by removing the trays and the rack. 3 Remove the filters at the bottom of the chamber (bottom/center) by unscrewing them. 4 Rinse the filters with tap water. 5 Insert the filters in their original position, screwing them into place.</p>	<p>LEXA DOOR GASKET</p> 			 <p>W&amp;H Lexa PM Kit</p>
1200 Cycle or Yearly	<p><b>Replace the Door Gasket, Chamber Filters, &amp; Hepa Filter</b></p> <p><i>Notice: The door gasket needs to be replaced every 1200 sterilization cycles or once a year, whichever come first. A replacement message alerts when replacement is due. If the consumables is replaced prior to the message, you have to reset the consumable cycle counter.</i></p> <p>1 Open the chamber door. 2 Remove the used door gasket by hand. 3 Carefully clean the seal seat and the inside face of the chamber door.</p>	<p>LEXA PM KIT</p>			 <p>W&amp;H Lexa PM Kit</p>



W&H Lexa

StatimUSA.com - Free Technical Support Line 704-966-1650 Option 3



FREE TECH SUPPORT



W&H Lexa User Manual