


# VistaPure Filters

## Original Design & Updated Design

**The Original Design**



**The Updated Design**



Why the new look?

We've made the new labels easier to read and reoriented the design so it's simple to see at a glance that the flow direction is correct when installing a new filter. Plus, the newly numbered, color-coded filter positions make replacement as easy as 1-2-3-4-5!

These new filters (shipped with all VistaPure systems manufacture after April 1, 2016) are compatible with older VistaPure systems.

Item #	Part #	Description	Filter Service Life	Old Part # Before 4/1/16
1	R3251	Prefilter	1 year max	R5633
2	R3252	Deionization Filter	4 -18 months	R5662
3	R3253	Deionization Filter	4 -18 months	R5662
4	R3254	Hyperfiltration Element	3 – 5 years	R3080
5	R3255	Polishing Filter	1 year max	R5633

**Note: These filters may be used on all VistaPure units**